# Dräger Alcotest<sup>®</sup> 6810 med

Atemalkohol-Messgerät  
Breath Alcohol Monitor

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# Contents

<b>For Your Safety and That of Your Patients</b> .....	4
<b>Intended use</b> .....	6
<b>What's what?</b> .....	7
Case .....	7
Measuring instrument .....	8
<b>Preparation</b> .....	9
Inserting the batteries .....	9
Fixing the hand strap .....	11
Symbols and readings .....	12
Fitting the mouthpiece .....	14
Switching on and off .....	15
Switch on .....	15
Switch off .....	15
Automatic switch off .....	15
<b>Performing the Measurement</b> .....	16
Conditions .....	16
Requirements of the person to be tested .....	16
Measuring .....	17
Automatic measurement .....	18
Error in delivering the breath sample .....	19
Removing the mouthpiece .....	20
Manual measurement .....	21
<b>Functional Test</b> .....	22
<b>Menu</b> .....	23
Activation .....	23
Operation .....	23

Symbole in der Anzeige .....	24	Symbols on the display .....	24
Menü .....	24	Menu .....	24
Funktion .....	24	Function .....	24
Menüpunkte .....	25	Menu items .....	25
<b>Fehler .....</b>	<b>26</b>	<b>Error .....</b>	<b>26</b>
<b>Technisches Handbuch .....</b>	<b>26</b>	<b>Technical Handbook .....</b>	<b>26</b>
<b>Instandhaltung .....</b>	<b>27</b>	<b>Maintenance .....</b>	<b>27</b>
Instandhaltungsintervalle .....	27	Maintenance intervals .....	27
Lagerung .....	27	Storage .....	27
Reinigen / Desinfizieren .....	28	Cleaning / Disinfecting .....	28
Batterien auswechseln .....	28	Changing batteries .....	28
Gerät entsorgen .....	30	Disposal of the instrument .....	30
<b>Störung, Ursache, Abhilfe .....</b>	<b>31</b>	<b>Fault, Cause, Remedy .....</b>	<b>32</b>
<b>Technische Daten .....</b>	<b>33</b>	<b>Technical data .....</b>	<b>33</b>
<b>Bestell-Liste .....</b>	<b>36</b>	<b>Order list .....</b>	<b>36</b>
Zubehör .....	36	Accessories .....	36
Ersatzteile .....	37	Spare parts .....	37
<b>Stichwortverzeichnis .....</b>	<b>38</b>	<b>Index .....</b>	<b>38</b>

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# For Your Safety and That of Your Patients

## **Strictly follow the instructions for use**

Any use of the instrument requires full understanding and strict observation of these instructions.  
The instrument is only to be used for the purposes specified herein.

## **Maintenance**

The instrument must be inspected and serviced regularly by trained service personnel at six monthly intervals (and a record kept).  
Repair of the instrument may only be carried out by trained service personnel.  
We recommend that a service contract be obtained with DrägerService and that all repairs also be carried out by them.  
Only authentic Dräger spare parts may be used for maintenance.  
Observe chapter "Maintenance Intervals".

## **Accessories**

Only the accessories specified in the order list may be used.

## **Not for use in areas of explosion hazard**

Alcotest<sup>®</sup> 6810 med is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

## **Safe connection with other electrical equipment**

Electrical connections to equipment which is not listed in these Instructions for Use should only be made following consultations with the respective manufacturers or an expert.

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Dräger Safety AG & Co. KGaA

### **Liability for proper function or damage**

The liability for the proper function of the instrument is irrevocably transferred to the owner or operator to the extent that the instrument is serviced or repaired by personnel not employed or authorized by Dräger Safety or if the instrument is used in a manner not conforming to its intended use.

Dräger Safety cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Dräger Safety are likewise not modified by the recommendations given above.

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## Intended use

For rapid measurement of breath alcohol concentration – for example for diagnosis of vascular injuries in the urology during transurethral resections of the prostate, TUR-P, or in the gynaecology during transcervical resections of the endometrium, TCRE, after previous application of an ethanol containing rinsing solution, at the emergency reception area for testing patients for previous alcohol consumption, or for the detection of alcohol consumption before methadone substitution.

**Do not use mobile phones within a distance of 10 metres from the instrument.**

**Mobile phones and portable or mobile HF communication equipments can cause interference to electrical and electronic medical appliances, thereby putting patients at risk.\***

**Do not use the instrument together with flammable gases or flammable anaesthetic agents. Fire hazard.**

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
---

® Alcotest is a registered trademark of Dräger.

\* Dräger medical devices comply with the interference immunity requirements of the specific standards for the products. However, depending on the design of the mobile phone and situation of use, field strengths may occur in the immediate environment of a mobile phone that exceed the values of the standards quoted and therefore may cause interference.

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

# Preparation

## Inserting the batteries

Before initial operation or if the battery symbol on the display indicates "empty".

### Only for NiMH-batteries

NiMH-batteries can be charged with a special power supply inside of the Alcotest 6810 med.

- Connect the special power supply (see page 36) with the charging/communication port.
- The charging process is indicated by a pulsing display of the charging symbol.


**Caution: Only use the power supply listed in the "Accessories" on page 36. Only charge NiMH-batteries.**

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00423855\_1.epa

## Batteries

- Gently press onto the lock of the battery compartment and remove cover.
- Observe correct polarity when inserting the batteries!
- Close battery cover.

When using NiMH batteries:


- Close special battery cover (see "Order list" on page 36) and fix the screw.

## Erforderliche Batterien

2 Stück 1,5 V Typ Alkaline (Mignon, LR6, AA)

oder

Akku 2 Stück – 1,2 V NiMH  
(Mignon, LR6, AA)

00423855\_1.epa

## Batteries to be used


2 pieces 1.5 V Type Alkaline (Mignon, LR6, AA)

or

rechargeable battery 2 pieces – 1.2 V NiMH  
(Mignon, LR6, AA)

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## Fixing the hand strap

- Insert the smaller loop into the hole in the housing.
- Thread the longer loop through the eyelet and pull through.

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



## Symbols and readings

The measured breath alcohol value is given in mg/L (= milligram ethanol per liter breath air).

### Symbols

In these instructions for use several symbols are used corresponding to the messages indicated on the display:

### Battery symbol

Indicates the battery charging status:

Full battery

Battery capacity 2/3


Battery capacity 1/3

Empty battery

### Note:

The symbols on the type plate and under the battery flap mean: **Attention! Please consult the user instructions!**  
**Before changing batteries, please read the instructions.**

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## Buttons


"OK"-button  
Confirms entries or selected functions.  
Button to switch instrument on and off.


"DOWN" and "UP"  
Setting of values and selection of menu functions.

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
006238465\_1.epa

## Fitting the mouthpiece

- Open packaging. To keep the mouthpiece hygienic, retain the portion to come in contact with mouth covered by packaging until mouthpiece is correctly fitted to the instrument.
- With the grooved markings facing up, apply the mouthpiece diagonally, shift it upwards and press it into mouthpiece receptacle until audible sound confirms correct engaging.
- Remove packaging completely.

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007238465\_1.epa

Mouthpiece can be used with orientation to right or left side.

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
## Switching on and off

### Switch on

- Keep  - button pressed for approx. 1 second until start window is displayed. A selftest is carried out. Device settings are displayed.

If the date for next calibration is expired:

- Unit has to be calibrated (see "Fault, Cause, Remedy" on page 32).
  - The use for measuring is still possible.

If the temperature inside the instrument is outside of the permissible temperature range, a measurement is not possible:

- Warm up, respectively cool down the instrument.

### Switch off

- Keep  - button pressed.


The instrument switches off after 2 seconds.

### Automatic switch off

The instrument switches off automatically 4 minutes after the last use.

When the batteries are discharged the instrument switches off automatically – the battery symbol is flashing.

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# Performing the Measurement

## Conditions

- A new mouthpiece is to be used for every person tested.
- No alcohol and solvent vapours (disinfectants) in the ambient air.
- Do not apply too high alcohol concentrations to sensor, e.g. by rinsing the mouth with high-grade alcohols shortly before the measurement. By doing so the sensor lifetime is reduced.
- Keep at a distance of at least 10 meters from the antennas of mobile phones and transmitting stations.

## Requirements of the person to be tested

- Maintain a waiting period of at least 15 minutes after drinking alcohol!  
Residual alcohol in the mouth can distort the measurement. Such residues may be left by aromatic drinks (e.g. fruit juice), mouth sprays containing alcohol, medicines and drops, as well as by burping and vomiting. In this case, wait for 15 minutes, too. Rinsing the mouth with water or non-alcoholic drinks does not reduce this waiting period!
- Before sampling the person to be tested must breathe normally. Quick inhaling and exhaling through the mouth must be avoided.

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After 2 seconds "**WAIT**" is indicated on the display.

After approximately 6 seconds "**READY**" is displayed and a short audible signal is given. At the same time the current test number is indicated. A green light flashes.

The instrument is ready for measurement.

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
01923895\_L\_de.npa

## Automatic measurement

During the automatic measurement the sample is automatically taken by the instrument when the person being tested has supplied a sufficient breath sample into the mouthpiece and the end of the breath sample has been detected by the instrument.

- Ask the test subject to blow evenly into the mouthpiece as long as possible without interruption.

During delivery of the breath sample a continuous tone sounds, the green light is lit, and the display shows: "**Blow**".

After delivery of the breath sample the green light will cease, the continuous tone will stop and a click will be heard.

The display shows: "**WAIT Analyzing**"


### Measurement result

After 5 to 25 seconds (depending on the instrument temperature and the measured concentration) the measured value is shown on the display.

The value is shown with the respective measuring unit (mg/L = milligram ethanol per liter breath air).

Documentation of the result:

While the result is shown on the display, a print-out can be produced with an appropriate Dräger printer.

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
## Further measurement

- By pressing the  - button a new measurement is prepared. During the recovery period of the sensor (see "Technical data" on page 33) the display shows: "WAIT"

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
### Note:

**Switching the instrument off in the meantime does not reduce the waiting period!**

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
If the sampling volume is too low the display shows: "**Insufficient Volume.**" The horn sounds briefly and a red light flashes.

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To repeat press the  - button. The instrument is ready for measurement after approx. 4 seconds.

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
When exhaling improperly, e.g. by abruptly stopping to exhale or sucking at the end of the sampling, the display shows: "Blow Interruption" The horn sounds briefly and a red light flashes.

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
Press the  - button to repeat the measurement. Time period until the instrument is ready for measurement again: 4 to 30 seconds.

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## Removing the mouthpiece

- Push the mouthpiece out of the mouthpiece holder.
- Dispose of the mouthpiece in accordance with local regulations.

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## Manual measurement

During the manual measurement the sampling system is manually activated by the user of the instrument.

This enables the use of the instrument also with persons being tested who are limited in their cooperation capabilities and are not able to blow into the mouthpiece on request or to deliver a sufficient breath sample.

- Hold the **instrument without a mouth-piece** into the air to be analyzed.
- Press the  - button shortly.

When the result is displayed a measured alcohol concentration is indicated additionally.

- Press the  - button for a new measurement.

NO  
ALKOHOL  
0.00 mg/L

02523815\_1\_01.de.qepa

### Note:

For the interpretation of the measuring result of a manual measurement the following has to be considered.

For a very small delivered breath volume the analysed breath air is a mixture of lung air and air from the upper respiratory tract. Therefore, the measured concentration is smaller compared to the blood alcohol concentration.


For a manual measurement **without** mouthpiece additionally a mixture of breath air with ambient air can occur, whereby the breath air is diluted and therefore the measured concentration is smaller than the concentration in the breath air.

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
# Functional Test

- Every week perform a functional test.
- Place the Alcotest® 6810 med ready for measurement without mouthpiece onto the opening of a bottle containing a solution with a concentration of approx. 1 g ethanol per 1 litre water (see Order List – Accessories, page 36) or approx. 1% ethanol rinsing solution.  
Alternatively drop some drops of the solution on an absorbent paper and place the instrument without mouthpiece above the paper.
- Press  - button shortly.
- The instrument performs a measurement and displays the result.  
As soon as the sensitivity of the instrument for alcohol has been verified the display shows "ALCOHOL", the red light is lit.
- Press  - button shortly; the instrument is ready for use.

If the display shows "NO ALCOHOL" and the green light is lit:

- Repeat functional test.
- If "NO ALCOHOL" shows again:  
Calibrate the instrument; by trained service personnel for example by the service of Dräger Safety.

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# Menu

## Activation

Switch instrument on. When "WAIT" or "READY" appears on the display the menu functions can be called up.

- For menu press button  shortly.

## Operation

- Select functions using  and  buttons.
- The selected function is highlighted.
- Activate the selected function by pressing  - button.
  
- Selection "  **return to measurement**" and pressing the  - button brings the instrument back into the measuring mode.

If no button is pressed within 120 seconds, the instrument automatically returns to the measuring mode.

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### Last tests

Alcotest<sup>®</sup> 6810 med stores the last 250 test results.

Selection of the tests with  and  . Back with  - button.

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### Last cal. date

Display of last calibration date. Back with  - button.


### Date & time

Display of present date and time. Back with  - button.

### Set date

Display and adjustment of date and date format. Change by pressing  and  , acknowledge by pressing  - button.

### Set time


Display and adjustment of time, time format, and automatic change between summer and winter time. Change by pressing  and  , acknowledge by pressing  - button.

### Set language

Selection of display language.

### Set contrast


Change by pressing  and  . Back with  - button.

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## Error

In case of internal error, the following error message is displayed (see page 32).

If the failure still occurs after the instrument is switched on and off several times resp. after removal of the batteries (at least 2 seconds) and re-inserting, contact Dräger Safety.

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## Technical Handbook

Additional information, for example for calibration or other options of measuring unit, is contained in the Technical Handbook (see order List – Accessories, page 36).

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- Every week, page 22.

### Safety inspection:

- Every six months, including calibration by trained service personnel, e. g. by Dräger Safety, including test of calibration.

## Storage

- The Dräger Alcotest<sup>®</sup> 6810 med consumes a small amount of current even when switched off. Therefore remove the batteries from the instrument during long periods without use.

### Note: In case of long storage periods

Remove batteries from instrument, if the instrument is not used for more than 6 months to avoid damage to the instrument caused by leaking battery acid.

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## Cleaning / Desinfecting


- Do not submerge in liquid – liquid must not enter the sockets.
- Do not sterilize in hot steam.
- Do not sterilize with ethylene oxide.

To clean:

- wipe with a disposable cloth.

To disinfect:

- wipe with a disposable cloth soaked in isopropyl alcohol.
- Wait at least 15 minutes until next measurement.

00423895\_1.epa

## Changing batteries

If the battery symbol on the display indicates "empty", change the batteries, resp. charge the rechargeable batteries (page 9).

- Gently press onto the lock of the battery compartment and remove cover.

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
00523856\_1.epa

- Remove used batteries and insert new ones or externally charged rechargeable batteries – consider correct polarity.
- Close battery cover.


#### **Required batteries:**

2 pieces 1.5 V Type Alkaline (Mignon, LR6, AA)

or

rechargeable battery 2 pieces 1.2 V NiMH  
(Mignon, LR6, AA)

For continuity of use of the instrument after discharge of the batteries it is recommended always to have replacement batteries available.

00523856\_1.epa

#### **Caution:**

##### **Used batteries**

- Do not expose to fire,
- Do not force open. Danger of an explosion!
- Do not recharge alkaline batteries!


Only dispose of as special waste, in accordance with local waste disposal regulations.

Information can be obtained from environmental agencies and from waste disposal companies.

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08223885\_1.qxp

## Disposal of the instrument

EU-wide regulations for the disposal of electric and electronic appliances which have been defined in the EU Directive 2002/96/EC and in national laws are effective from August 2005 and apply to this device.

Common household appliances can be disposed of using special collecting and recycling facilities. However, as this device has not been registered for household usage, it must not be disposed of through these means. The device can be returned to your national Dräger Safety Sales Organization for disposal. Please do not hesitate to contact the above if you have any further questions on this issue.

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## Fault, Cause, Remedy

Fault	Cause	Remedy
Instrument switches off automatically.	Instrument has been left ready for measurement for more than 4 minutes.	Switch on again.
	Batteries discharged.	Replace batteries, page 28, resp. charge the rechargeable batteries.
Instrument leaves menu-mode automatically and shows measuring mode.	Instrument has an automatically switch off of the menu after 120 sec. without use.	Reselect menu item.
Instrument cannot switched on.	Batteries discharged.	Replace batteries, resp. charge the rechargeable batteries.
Insufficient volume, blowing interruption.	Testperson not blowing hard or evenly enough.	Ask to blow harder and more evenly.
Service date expired.	The service period has expired.	Perform safety inspection, by trained service personnel, e. g. from Dräger Safety.
No text in the display.	Instrument error, eventually caused by electrostatic discharge.	Remove batteries (at least 2 sec.) and insert them again, if the fault reappears call Service of Dräger Safety.

Fault in display	Cause	Remedy
"Accuracy check error"	Accuracy check not successful.	Repeat accuracy check, if fault reappears, call Service of Dräger Safety.
	Measured concentration exceeds measuring limits.	Use lower concentrations.
"Internal System ERROR" with service code	Instrument error, no measuring function possible.	Remove batteries (at least 2 sec.) and insert them again, if the fault reappears call Service of Dräger Safety.
"Service date expired"	The service period has expired.	Perform safety inspection, by trained service personnel, e. g. from Dräger Safety.
"Battery charging error"	No rechargeable NiMH batteries in the instrument.	Replace batteries with rechargeable NiMH batteries. Start charging.

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## Technical data

<b>Measuring principle</b>	Electrochemical sensor
<b>Measurement range</b>	
Breath alcohol concentration	0.00 bis 2.50 mg/L (Ethanol-Mass per breath volume at 34 °C and 1013 hPa.)

### Display range

Display range is identical to measuring range.  
If the result is exceeding the measuring range the display will show three arrows "↑ ↑ ↑" as overflow indication.

### Ambient conditions

For operation	-5 to 50 °C (23 to 122 °F) 10 to 100 % r.h. (non condensing) 600 to 1400 hPa
For storage	-40 to 70 °C (-40 to 158 °F)

**Sample conditions**  
(automatic measurement)

Breath volume	> 0.3 L
Breath flow	>3 L/min

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**Measurement accuracy**

Standard deviation of reproducibility with an ethanol standard

up to 0.50 mg/L	0.008 mg/L
≥ 0.50 mg/L	1.7 % of measured value

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keine

**Cross sensitivity to anesthetic gases**

none

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1. Messung	ca. 6 s
------------	---------

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0,25 mg/L	20 s
0,50 mg/L	40 s
1,00 mg/L	80 s
> 1,5 mg/L	120 s

**Intervals before readiness for operation**

1st measurement	approx. 6 s
-----------------	-------------

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## Order list

Name and Description	Part no.
<b>Dräger Alcotest® 6810 med</b> (instrument, 3 mouthpieces, batteries, hand strap and case)	<b>83 19 080</b>
<b>Accessories</b>	
<b>Set rechargeable battery</b> , for use with rechargeable NiMH-batteries (special battery cover, 2 rechargeable NiMH batteries)	<b>83 18 565</b>
<b>Case for Dräger Alcotest® 6XXX</b>	<b>83 17 920</b>
<b>Leather pouch for Dräger Alcotest® 6XXX (black)</b> (belt fixation)	<b>83 17 911</b>
<b>Leather pouch for Dräger Alcotest® 6XXX (yellow)</b> (Instrument can be operated while staying in the pouch) with pocket for 4 mouthpieces, no belt fixation	<b>83 17 931</b>
<b>"AlcoView"</b> PC-programm for communication with the Alcotest® 6810	<b>83 17 915</b>
<b>PC-connection cable</b>	<b>83 17 593</b>
<b>Power supply</b> (current source) for charging NiMH batteries inside the Alcotest® 6810	<b>83 16 991</b>
Technical Handbook	<b>90 23 833</b>
<b>Accessories for functional test:</b>	<b>67 28 838</b>
Ethanol standard solution 1.21 g/L (1 x 500 mL)	

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<b>Name</b> and Description	<b>Part no.</b>
<b>Spare parts</b>	
<b>Mouthpieces</b> with non-return valve (100 pieces)	<b>68 11 055</b>
<b>Mouthpieces</b> with non-return valve (250 pieces)	<b>68 11 060</b>
<b>Mouthpieces</b> with non-return valve (1000 pieces)	<b>68 11 065</b>
<b>Alkaline battery</b> (1 piece)	<b>13 35 804</b>
<b>Rechargeable battery NiMH</b> (1 piece)	<b>83 17 912</b>

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Alkaline battery supply	35
Ambient conditions	33
Breath alcohol concentration	33
Breath volume	33
Calibration date	25
Changing Batteries	28
Classification	35
Cleaning	28
Contrast	25
Date & time	25
Dimensions	35
Display range	33
Documentation	18
Electrical supply	35
Functional test	22, 34
Intervals before readiness for operation	34
Language selection	25
Last tests	25
Maintenance Intervals	27
Measurement accuracy	34
Measurement range	33
Measuring principle	33
Order List	36

Probenahme .....	33	Reproducibility .....	34
Reinigen .....	28	Safety inspection .....	27
Sicherheitstechnische Kontrolle .....	27	Sample conditions .....	33
Sondermüll .....	29	Sensitivity drift .....	34
Sprache wählen .....	25	Set date .....	25
Stromversorgung .....	35	Set time .....	25
Technische Daten .....	33	Special waste .....	29
Technisches Handbuch .....	26	Storage .....	27
Uhrzeit .....	25	Technical Data .....	33
UMDNS-Code .....	35	Technical Handbook .....	26
Umgebungsbedingungen .....	33	UMDNS-Code .....	35
Verbrauchte Alkali Batterien .....	29	Used alkaline batteries .....	29
Wartezeit für Messbereitschaft .....	34	Weight .....	35

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Directive 93/42/EEC  
concerning Medical Devices

**Dräger Safety AG & Co. KGaA**

Revalstraße 1

D-23560 Lübeck

Tel. (0451) 8 82- 0

Fax (0451) 8 82- 20 80

[www.draeger.com](http://www.draeger.com)

**Dräger Safety AG & Co. KGaA**

Revalstraße 1

D-23560 Lübeck

Germany

Tel. +49 451 8 82- 0

Fax +49 451 8 82- 20 80

[www.draeger.com](http://www.draeger.com)

90 23 855 - GA 4754.551 de/en

© Dräger Safety AG & Co. KGaA

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90 23 855 - GA 4754.551 de/en

© Dräger Safety AG & Co. KGaA

3rd edition - January 2007

Subject to alteration